


FITTING INSTRUCTIONS:

1

**Step 1**

Hold respirator in your hand as shown, with the nosepiece at your fingertips and the headbands hanging freely below your hand.

2

**Step 2**


Place respirator under your chin with the nosepiece up. Pull the top headband to top back of head.

3

**Step 3**

Pull the bottom headband overhead and place around neck below ears.

4

**Step 4**

Adjust headbands for comfortable fit. Push the nosepiece inward with your fingertips to conform to the shape of your nose.

5

**Step 5**

Cup both hands over the mask, checking for tightness, Inhale sharply, negative pressure should be felt in the mask, Re-adjust as necessary.

Caution:

The respirator should be fit checked for each wearing. Place both hands completely over the respirator and exhale. If air leaks around nose or the edges, adjust the nosepiece and/or headbands until a good fit is achieved.